

Lot-No.

MANUAL - one step

Ref. K102

Expiry date: 1 year

100 Tests (Ready to use PCR kit)

STORE AT -20°C

Japanese Encephalitis Virus (JEV) – one step PCR

-Only for in vitro use-
-Only for veterinary use-
-Humans, only for research use-
-To be used by a technical person-

Principle and use

This amplification kit has been manufactured by *Genekam Biotechnology AG*, Germany to detect *Japanese Encephalitis virus (JEV)* in two steps: 1. Flavivirus in one step 2. JEV in the other step.

This kit needs RNA which can be isolated from blood, serum, tissue, mosquitoes, and any body fluid. Kindly use good methods to isolate the RNA. Kindly take common safety laboratory precautions during working. ***Please use gloves during work. Proceed clean and carefully otherwise you may cause contamination problems. Do not touch other objects like pens, chairs etc. during Part 1.***

IMPORTANT: we added cotton or sponge in the lid of container of the kit, to avoid damage during transportation. Please remove this cotton or sponge from the lid of each container before storage.

Composition:

It contains the following (**WARNING! THAW THE TUBES SLOWLY: NEVER THAW IN HEATING BLOCK OR WITH HEAT FROM HAND**):

- Tube A (for Flavivirus) (2 tubes)
- Tube B (3 tubes)
- Tube C (for JEV) (2 tubes) : Sometime this tube is also called H tube.
- Tube Y (1 tube)
- positive (+ve) control (D1) (1 tube): **it should be stored at -20°C. It is for yellow fever.**
- positive (+ve) control (D3) (1 tube): **it should be stored at -20°C. It is for JEV.**
- Negative (-Ve) Control (tube D2) (1 tube)
- Marker (tube E) (1 tube) (max 1000bp) : 100, 150, 200, 300, 400, 500, 600, 700, 800, 900, 1000bp
- Dye (tube F) (1 tube)

Please check them before you start.

Equipment needed:

- PCR thermocycler
- Laboratory centrifuge
- UV platform
- UV safety goggles
- microtubes (0.2ml)
- Pipette-tips with and without filter (20µl, 5µl & 1µl)
- Pipettes (quality pipettes)
- Gel Agarose chamber
- Power supply
- Paper
- Pen
- Agarose (good quality)
- Staining (Ethium Bromide)
- TAE buffer 1x

- Ice
- Vortexer



Procedure:

PART 1 – it is a one step PCR.


STEP A


1. Kindly thaw **one tube** each: Y, A, B, C, D1, D2, D3, E and F (don't thaw in hot block). After thawing, kindly put the tubes at 4°C. If the kit is not in use, store them at -20°C.

2. Mark your microtubes with a sample number, +ve Control and -ve Control. Please divide tubes in two sets. One set is marked with word (Flavi), where the other set is marked word (JEV).





3. Add 1µl of tube Y to each tube of both sets.






4. Add 10µl of B to each micro tube of both sets. Avoid to touch the wall of the micro tubes.




5. Add 7µl of A to each tube of the set Flavi (avoid to touch the wall of the microtubes).




TIP: to save time and money, you can calculate how much chemicals you need to run the test. You want to run 10 test, i.e. you need $10 \times 7\mu\text{l}$ of A = $70\mu\text{l}$ of A + $10 \times 10\mu\text{l}$ of B = $100\mu\text{l}$ of B + $10 \times 1\mu\text{l}$ of Y = $10\mu\text{l}$ of Y = $180\mu\text{l}$ in total. From this, $18\mu\text{l}$ can be added to each microtube.

6. Add $7\mu\text{l}$ of C to each tube of the set JEV (avoid to touch the wall of the microtubes).




TIP: to save time and money, you can calculate how much chemicals you need to run the test. You want to run 10 test, i.e. you need $10 \times 7\mu\text{l}$ of C = $70\mu\text{l}$ of C + $10 \times 10\mu\text{l}$ of B = $100\mu\text{l}$ of B + $10 \times 1\mu\text{l}$ of Y = $10\mu\text{l}$ of Y = $180\mu\text{l}$ in total. From this, $18\mu\text{l}$ can be added to each microtube.


7. Add $2\mu\text{l}$ of your RNA with pipette tip with filter to each micro tube of both sets according to your label except +Ve and -Ve (Avoid touching the wall). Use every time a new pipette tip (for each sample)! Mix it thoroughly.




8. Use new pipette tip with filter. Add $2\mu\text{l}$ of D1 to positive control tube (Flavi). Use a new pipette tip. Mix it.




8.1 Use new pipette tip with filter. Add 2µl of D3 to positive control tube (JEV). Use a new pipette tip. Mix it.



9. Use a new pipette tip. Add 2µl of -Ve (tube D2) to -Ve Control (Flavi). Mix it.





9.1 Use a new pipette tip. Add 2µl of -Ve (tube D2) to -Ve Control (JEV). Mix it.



10. Centrifuge all tubes for 20 sec. for 8000 rpm (this is not necessary but it is better).

11. Run the program of your thermocycler as followings:
Kindly check whether you have added everything correctly as the level of the volume of each micro tube must be almost the same.



Now program your PCR machine as follows.

1. 3600 seconds at 42°C
600 seconds at 70°C
 2. 300 seconds at 37°C
 3. A. 60 seconds at 94°C
B. 60 seconds at 53°C
C. 60 seconds at 72°C
- } 30 cycles


Before you start the PCR program, kindly check whether tubes are closed properly. **Microtubes must be in contact with metal block** (very important!). There should be no air or loose contact with metal block of thermocycler.

12. After step 11 is finished take out the microtubes and centrifuge for a while.


Now go to STEP B.

STEP B

1. Prepare the gel Agarose 2% in TAE (1x) buffer.
2. Let the Gel dry and add this TAE (1x) buffer in gel chamber.
3. Take the tube E (marker). Make ready to use for gel electrophoresis.
4. After the PCR step is finished, now you can prepare for gel Agarose electrophoresis. Take 2µl of Dye (tube F) and add to each micro tube (with the same number as your PCR microtubes including +Ve, -Ve controls).



5. Add 10µl of marker (tube E: 100bp) to first and last lane of electrophoresis. (Kindly make lane plan on paper according to your probes in order to identify later and see the results).



6. Add 10µl of mix of step 4 to each lane of gel Agarose (between first and last lane). Change the pipette tip for each lane.

7. Run the gel for **55 min.** at **120 Volt**. It may vary.

8. Make staining solution ready.

9. Put the gel for 5-30 minutes staining solution (0.5µg/ml).

10. View the gel under UV light. UV light is dangerous for your eyes. Use UV goggles.

11. You must find the bands in positive control and no band in negative control. The interpretation of the results are available on extra page.

Tip: Gene sequencing is highly recommended to reconfirm the positive results.

If you should find any mistakes, please let us know. Thank you.

Suggestion:

This manual has been written specifically for beginners, hence persons with experience in PCR must use their experience to keep each step as small as possible e.g. you should calculate the amount of the needed chemicals, before starting with testing.

Last update: 31-01-2012
v1.2

Genekam Biotechnology AG

Dammstr. 31-33

47119 Duisburg

Germany

Tel. (+49) 203 / 555858-31,-32,-33

Fax (+49) 203 / 35 82 99

anfrage@genekam.de

<http://www.genekam.de>

K102: Japanese Encephalitis virus - RESULTS

1. In PCR product of Part A, you will see the following bands:

675 band =====> Yellow fever (control cDNA): *this is included as control in our kit.*

673 and 751 band => Japanese Encephalitis

760 and 834 band => West Nile virus

541 and 628 band => D2 (Dengue)

552 and 607 band => D3 (Dengue)

471 and 560 band => D4 (Dengue)

2. In PCR product of Part B, you will have a band of **142 bp**. This is specific for JEV.