

Lot.
Ref. SB0096

MANUAL

Expiry date: 1 year

STORE AT -20°C

dNTP 10mM each Total 500 ul (1 Tubes)

-Only for in vitro use-
-Only for research use-
-To be used by a technical person-

This contains 1 tubes. Each tube contains 500ul dNTP 10nM each. dNTP mix can be used in different applications like PCR master, RNA to cDNA conversion.

IMPORTANT: we added cotton or sponge in the lid of container of the kit, to avoid damage during transportation. Please remove this cotton or sponge from the lid of each container before storage.

Material safety details:

- there are no material safety details.
- Please use gloves during working.

QUALITY CONTROL:

- passed through quality control